



**WRIGHT
&
FILIPPIS**

BUTLER WHEELCHAIR LIFTS

Butler Wheelchair Lifts

**PROVIDE
SAFE,
CONVENIENT
ACCESS
FOR ALL**

From single family private residences to public commercial facilities, Butler Wheelchair Lifts provide safe, convenient access for everyone. Adaptable to nearly all installation applications, Butler Lifts are long-lasting, affordable and surprisingly easy to install.



As Easy To Install As It Is To Use

Butler Wheelchair Lifts are almost as easy to install as track lighting for Wright & Filippis. We will have your lift installed and running in a short period of time.

If door frames protrude, the lift's guide unit can be moved to clear it. The power unit can be located just about anywhere it's convenient and out of the way, such as through a wall.



Heavy gauge rails use just 4" of stairway on either side. Rail attaches directly to landing at top and bottom.

The Butler Lift Track

The lift attaches directly to the top and bottom stair treads - no anchoring to the wall is required. With a special built-in lowering system and power outage safety features, Butler Lifts provide maximum rider safety.

Strength Plus Beauty

Just because a wheelchair lift is a necessity doesn't mean it necessarily has to be an eyesore.

For example, all (and we do mean *all*) of the lift's working mechanisms have been cleverly concealed to provide a more attractive, sleeker appearance.

Offering a lifting capacity of 500 lbs., these lifts can be equipped with an end ramp and guard which is automatically retractable and locks firmly in place when the lift is in operation or at rest at the top of the stairs. And a built-in lowering system makes it easy for even a young child to clear the stairway in an emergency.

There's also no need to worry about power outages leaving the user stranded between floors. With a simple push of a lever, the lift will safely and gently move to the lower level in spite of the power failure.

Three sets of directional control switches are located within easy access of the user, and require constant pressure for additional safety.



Plus, Butler Wheelchair Lifts take up only 8" of stairway space when not in use. This unique narrow profile means increased safety for both the user and other members of the household.

The Details

Track

Heavy gauge rectangular tubing encasing trolley rollers, slack and broken chain safety, speed safety and roller chain.

Roller Chain

Used for suspension with a chain on both sides of the stairs. Standard rated strength is 3700 lbs. on both sides under a load of 500 lbs. plus weight of platform.

Controls

Three directional control switches are located on the car platform, at the upper landing under one of the tracks and at the lower landing at the end of the machinery cover towards the stairway.

Speed Of Travel

22 feet per minute.

Current

110 volts, single phase, 50/60 cycle, standard household current. Power cord prepared for inserting into receptacle. All control, limit and safety switches are on a safe 24 volt system.



The end ramp locks firmly into place in the up position when the lift is in use.

Easy-to-reach hand controls make the Butler Lifts simple to operate.

Unique hand control with emergency brake release lever allows the unit to drift down safely in the event of a power failure.



WRIGHT & FILIPPIS

First to Serve, First to Care...Since 1944

1-800-482-0222 or 248-829-8214

www.FirstToServe.com

With over 20 locations throughout Michigan, please call for the location nearest you.